

# CLEC4G Antibody, HRP conjugated



PACO63128

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds mannose, N-acetylglucosamine (GlcNAc) and fucose, but not galactose, in a Ca(2+)-dependent manner, in vitro.

**Gene ID:**

CLEC4G

**Uniprot**

Q6UXB4

**Synonyms:**

C-type lectin domain family 4 member G (Liver and lymph node sinusoidal endothelial cell C-type lectin) (LSEctin), CLEC4G

**Immunogen:**

Recombinant Human C-type lectin domain family 4 member G protein (66-293AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A