

TNFAIP8L2 Antibody, FITC conjugated



PACO63141

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a negative regulator of innate and adaptive immunity by maintaining immune homeostasis. Negative regulator of Toll-like receptor and T-cell receptor function. Prevents hyperresponsiveness of the immune system and maintains immune homeostasis. Inhibits JUN/AP1 and NF-kappa-B activation. Promotes Fas-induced apoptosis.

Gene ID:

TNFAIP8L2

Uniprot

Q6P589

Synonyms:

Tumor necrosis factor alpha-induced protein 8-like protein 2 (TIPE2) (TNF alpha-induced protein 8-like protein 2) (TNFAIP8-like protein 2) (Inflammation factor protein 20), TNFAIP8L2

Immunogen:

Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 2 protein (1-184AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A