

DDX47 Antibody, HRP conjugated



PACO63212

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in apoptosis. May have a role in rRNA processing and mRNA splicing. Associates with pre-rRNA precursors.

Gene ID:

DDX47

Uniprot

Q9H0S4

Synonyms:

Probable ATP-dependent RNA helicase DDX47 (EC 3.6.4.13) (DEAD box protein 47), DDX47

Immunogen:

Recombinant Human Probable ATP-dependent RNA helicase DDX47 protein (225-460AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A