

DVL1P1 Antibody



PACO63287

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

May play a role in the signal transduction pathway mediated by multiple Wnt genes.

Gene ID:

DVL1P1

Uniprot

P54792

Synonyms:

Putative segment polarity protein dishevelled homolog DVL1P1 (DSH homolog 1-like) (Segment polarity protein dishevelled homolog DVL-1-like) (Dishevelled-1-like), DVL1P1, DVL DVL1 DVL1L1


Immunogen:

Peptide sequence from Human Putative segment polarity protein dishevelled homolog DVL1P1 protein (152-173AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: PC-3 whole cell lysate. All lanes: DVL1P1 antibody at 1:1000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 74, 73 kDa. Observed band size: 74 kDa.