

# FGFR2 Antibody, FITC conjugated



PACO63309

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FGFR2

**Uniprot**

P21802

**Synonyms:**

Fibroblast growth factor receptor 2 (FGFR-2) (EC 2.7.10.1) (K-sam) (KGFR) (Keratinocyte growth factor receptor) (CD antigen CD332), FGFR2, BEK KGFR KSAM

**Immunogen:**

Recombinant Human Fibroblast growth factor receptor 2 protein (577-821AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A