

# HSPA12A Antibody



PACO63323


---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ul	extracellular exosome.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	HSPA12A
<b>Source:</b>	<b>Uniprot</b>
Rabbit	O43301
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Heat shock 70 kDa protein 12A, HSPA12A, KIAA0417
<b>Applications:</b>	<b>Immunogen:</b>
ELISA, WB	Recombinant Human Heat shock 70 kDa protein 12A protein (2-231AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
ELISA:1:2000-1:10000, WB:1:500-1:2000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in: SH-SY5Y whole cell lysate, 293 whole cell lysate, HepG2 whole cell lysate. All lanes: HSPA12A antibody at 1:1000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 75 kDa. Observed band size: 75 kDa.