

PAX8 Antibody, HRP conjugated



PACO63368

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

Gene ID:

PAX8

Uniprot

Q06710

Synonyms:

Paired box protein Pax-8, PAX8

Immunogen:

Recombinant Human Paired box protein Pax-8 protein (237-450AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A