

DUSP13 Antibody, HRP conjugated



PACO63508

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Dual specificity phosphatase that dephosphorylates MAPK8/JNK and MAPK14/p38, but not MAPK1/ERK2, in vitro. Exhibits intrinsic phosphatase activity towards both phospho-seryl/threonyl and -tyrosyl residues, with similar specific activities in vitro.

Gene ID:

DUSP13

Uniprot

Q9UII6

Synonyms:

Dual specificity protein phosphatase 13 isoform B (DUSP13B) (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase SKRP4) (Testis- and skeletal-muscle-specific DSP), DUSP13, DUSP13B TMDP

Immunogen:

Recombinant Human Dual specificity protein phosphatase 13 isoform B protein (1-198AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A