

# SLC39A10 Antibody, HRP conjugated



PACO63520

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May act as a zinc-influx transporter.

**Gene ID:**

SLC39A10

**Uniprot**

Q9ULF5

**Synonyms:**

Zinc transporter ZIP10 (Solute carrier family 39 member 10) (Zrt- and Irt-like protein 10) (ZIP-10), SLC39A10, KIAA1265 ZIP10

**Immunogen:**

Recombinant Human Zinc transporter ZIP10 protein (41-317AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A