

SLC39A10 Antibody, Biotin conjugated



PACO63522

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as a zinc-influx transporter.

Gene ID:

SLC39A10

Uniprot

Q9ULF5

Synonyms:

Zinc transporter ZIP10 (Solute carrier family 39 member 10) (Zrt- and Irt-like protein 10) (ZIP-10), SLC39A10, KIAA1265 ZIP10

Immunogen:

Recombinant Human Zinc transporter ZIP10 protein (41-317AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A