

IVNS1ABP Antibody, Biotin conjugated



PACO63534

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in cell division and in the dynamic organization of the actin skeleton as a stabilizer of actin filaments by association with F-actin through Kelch repeats. Protects cells from cell death induced by actin destabilization; Protects neurons from dendritic spines and actin filaments damage induced by the actin-destabilizing cytochalasin B when overexpressed. Activates Erk signaling pathway when overexpressed in cultured cell lines. May be a component of the cellular splicing machinery with a role in pre-mRNA splicing; may mediate the inhibition of splicing by NS/influenza virus NS1A protein. Functions as modifier of the AHR/Aryl hydrocarbon receptor pathway increasing the concentration of AHR available to activate transcription.

Gene ID:

IVNS1ABP

Uniprot

Q9Y6Y0

Synonyms:

Influenza virus NS1A-binding protein (NS1-BP) (NS1-binding protein) (Aryl hydrocarbon receptor-associated protein 3), IVNS1ABP, ARA3 FLARA3 KIAA0850 NS1 NS1BP

Immunogen:

Recombinant Human Influenza virus NS1A-binding protein (401-642AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A