

CD4-like protein Antibody, Biotin conjugated



PACO63566

Product Information

Size:

50ul

Reactivity:

Grass carp

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CD4-like protein

Uniprot

C7FF20

Synonyms:

CD4-like protein

Immunogen:

Peptide sequence from Grass carp CD4-like protein (52-72AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A