

grxD Antibody



PACO63615

Product Information

Size:

50ul

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Monothiol glutaredoxin involved in the biogenesis of iron-sulfur clusters.

Gene ID:

grxD

Uniprot

P0AC69

Synonyms:

Glutaredoxin 4, Grx4, Monothiol glutaredoxin, grxD, ydhD


Immunogen:

Recombinant Escherichia coli Glutaredoxin 4 protein (24-115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: grxD antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 37 kDa. Observed band size: 35 kDa.