

COL6A1 Antibody, Biotin conjugated



PACO63646

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Collagen VI acts as a cell-binding protein.

Gene ID:

COL6A1

Uniprot

P12109

Synonyms:

Collagen alpha-1(VI) chain, COL6A1

Immunogen:

Peptide sequence from Human Collagen alpha-1(VI) chain protein (64-80AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A