

TMEM131L Antibody, FITC conjugated



PACO63665

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Isoform 1: Membrane-associated form that antagonizes canonical Wnt signaling by triggering lysosome-dependent degradation of Wnt-activated LRP6. Regulates thymocyte proliferation.

Gene ID:

TMEM131L

Uniprot

A2VDJ0

Synonyms:

Transmembrane protein 131-like, TMEM131L, KIAA0922

Immunogen:

Peptide sequence from Human Transmembrane protein 131-like protein (963-978AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A