

LGALS4 Antibody, HRP conjugated



PACO63672

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Galectin that binds lactose and a related range of sugars. May be involved in the assembly of adherens junctions.

Gene ID:

LGALS4

Uniprot

P56470

Synonyms:

Galectin-4 (Gal-4) (Antigen NY-CO-27) (L-36 lactose-binding protein) (L36LBP) (Lactose-binding lectin 4), LGALS4

Immunogen:

Recombinant Human Galectin-4 protein (1-323AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A