

TIAL1 Antibody, Biotin conjugated



PACO63718

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

RNA-binding protein. Possesses nucleolytic activity against cytotoxic lymphocyte target cells. May be involved in apoptosis.

Gene ID:

TIAL1

Uniprot

Q01085

Synonyms:

Nucleolysin TIAR (TIA-1-related protein), TIAL1

Immunogen:

Recombinant Human Nucleolysin TIAR protein (9-277AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A