

VIPR2 Antibody, Biotin conjugated



PACO63734

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is a receptor for VIP as well as PACAP-38 and -27, the activity of this receptor is mediated by G proteins which activate adenylyl cyclase. Can be coupled to phospholipase C.

Gene ID:

VIPR2

Uniprot

P41587

Synonyms:

Vasoactive intestinal polypeptide receptor 2 (VIP-R-2) (Helodermin-preferring VIP receptor) (Pituitary adenylate cyclase-activating polypeptide type III receptor) (PACAP type III receptor) (PACAP-R-3) (PACAP-R3) (VPAC2), VIPR2, VIP2R

Immunogen:

Recombinant Human Vasoactive intestinal polypeptide receptor 2 protein (24-126AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A