

CD24 Antibody, Biotin conjugated



PACO63738

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Modulates B-cell activation responses. Signaling could be triggered by the binding of a lectin-like ligand to the CD24 carbohydrates, and transduced by the release of second messengers derived from the GPI-anchor. Promotes AG-dependent proliferation of B-cells, and prevents their terminal differentiation into antibody-forming cells.

Gene ID:

CD24

Uniprot

P25063

Synonyms:

Signal transducer CD24 (Small cell lung carcinoma cluster 4 antigen) (CD antigen CD24), CD24, CD24A

Immunogen:

Peptide sequence from Human Signal transducer CD24 protein (28-46AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A