

PACO63751

Product Information

Size:

50ul

Reactivity:

Chicken anemia virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May act as a scaffold protein in virion assembly. May also play a role in intracellular signaling during viral replication.

Gene ID:

VP2

Uniprot

P69484

Synonyms:

Dual specificity protein phosphatase VP2 (EC 3.1.3.16) (EC 3.1.3.48), VP2


Immunogen:

Recombinant Chicken anemia virus Dual specificity protein phosphatase VP2 protein (1-216AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: VP2 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 62 kDa. Observed band size: 62 kDa.