

FABP9 Antibody



PACO63779

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

cytosol, triglyceride catabolic process.

Gene ID:

FABP9

Uniprot

Q0Z7S8

Synonyms:

Fatty acid, binding protein 9 (Testis lipid-binding protein) (TLBP) (Testis-type fatty acid, binding protein) (T-FABP), FABP9


Immunogen:

Recombinant Human Fatty acid, binding protein 9 protein (1-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of PC-3 cells with PACO63779 at 1:50, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).