

RAB39A Antibody, FITC conjugated



PACO63789

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in the maturation and acidification of phagosomes that engulf pathogens, such as S.aureus and M. tuberculosis. Plays a role in vesicular trafficking. Plays a role in the fusion of phagosomes with lysosomes. Negatively regulates LPS-induced autophagosome formation in macrophages possibly by implicating PI3K. May be involved in multiple neurite formation.

Gene ID:

RAB39A

Uniprot

Q14964

Synonyms:

Ras-related protein Rab-39A (Rab-39), RAB39A, RAB39

Immunogen:

Recombinant Human Ras-related protein Rab-39A protein (1-217AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A