

KCNV1 Antibody, Biotin conjugated



PACO63810

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potassium channel subunit that does not form functional channels by itself. Modulates KCNB1 and KCNB2 channel activity by shifting the threshold for inactivation to more negative values and by slowing the rate of inactivation. Can down-regulate the channel activity of KCNB1, KCNB2, KCNC4 and KCND1, possibly by trapping them in intracellular membranes.

Gene ID:

KCNV1

Uniprot

Q6PIU1

Synonyms:

Potassium voltage-gated channel subfamily V member 1 (Neuronal potassium channel alpha subunit HNKA) (Voltage-gated potassium channel subunit Kv8.1), KCNV1

Immunogen:

Recombinant Human Potassium voltage-gated channel subfamily V member 1 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A