

CACNG6 Antibody, Biotin conjugated



PACO63918

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thought to stabilize the calcium channel in an inactivated (closed) state.

Gene ID:

CACNG6

Uniprot

Q9BXT2

Synonyms:

Voltage-dependent calcium channel gamma-6 subunit (Neuronal voltage-gated calcium channel gamma-6 subunit), CACNG6

Immunogen:

Peptide sequence from Human Voltage-dependent calcium channel gamma-6 subunit protein (22-39AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A