

UTP11 Antibody, HRP conjugated



PACO63924

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in nucleolar processing of pre-18S ribosomal RNA.

Gene ID:

UTP11

Uniprot

Q9Y3A2

Synonyms:

Probable U3 small nucleolar RNA-associated protein 11 (U3 snoRNA-associated protein 11) (UTP11-like protein), UTP11, UTP11L

Immunogen:

Recombinant Human Probable U3 small nucleolar RNA-associated protein 11 protein (1-253AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A