

CHRND Antibody, HRP conjugated



PACO63952

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

Gene ID:

CHRND

Uniprot

Q07001

Synonyms:

Acetylcholine receptor subunit delta, CHRND, ACHRD

Immunogen:

Recombinant Human Acetylcholine receptor subunit delta protein (334-471AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A