

GIP Antibody, Biotin conjugated



PACO63958

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potent stimulator of insulin secretion and relatively poor inhibitor of gastric acid, secretion.

Gene ID:

GIP

Uniprot

P09681

Synonyms:

Gastric inhibitory polypeptide (GIP) (Glucose-dependent insulinotropic polypeptide) (Incretin hormone), GIP

Immunogen:

Peptide sequence from Human Gastric inhibitory polypeptide protein (113-131AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A