

PRAMEF6 Antibody, HRP conjugated



PACO63964

Product Information

Size:

50ul

Protein Background:**Reactivity:**

Human

Gene ID:

PRAMEF6

Source:

Rabbit

Uniprot

Q5VXH4

Isotype:

IgG

Synonyms:

PRAME family member 6, PRAMEF6

Applications:

ELISA

Immunogen:

Recombinant Human PRAME family member 6 protein (209-476AA).

Recommended dilutions:**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A