

OPCML Antibody, HRP conjugated



PACO63968

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds opioids in the presence of acid, c lipids; probably involved in cell contact.

Gene ID:

OPCML

Uniprot

Q14982

Synonyms:

Opioid-binding protein/cell adhesion molecule (OBCAM) (OPCML) (Opioid-binding cell adhesion molecule) (IgLON family member 1), OPCML, IGLON1 OBCAM

Immunogen:

Recombinant Human Opioid-binding protein/cell adhesion molecule protein (28-322AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A