

SEZ6L2 Antibody, HRP conjugated



PACO64016

Product Information

Size:

50ul

Reactivity:

Human

Source:

Seizure 6-like protein 2

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May contribute to specialized endoplasmic reticulum functions in neurons.

Gene ID:

SEZ6L2

Uniprot

Q6UXD5

Synonyms:

Seizure 6-like protein 2, SEZ6L2, PSK

Immunogen:

Recombinant Human Seizure 6-like protein 2 protein (28-372aa).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A