

Fgf2 Antibody, FITC conjugated



PACO64021

Product Information

Size:

50ul

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

Fgf2

Uniprot

P13109

Synonyms:

Fibroblast growth factor 2 (FGF-2) (Basic fibroblast growth factor) (bFGF) (Heparin-binding growth factor 2) (HBGF-2), Fgf2, Fgf-2

Immunogen:

Recombinant Rat Fibroblast growth factor 2 protein (10-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A