

ZNF346 Antibody



PACO64035

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Binds with low affinity to dsDNA and ssRNA, and with high affinity to dsRNA, with no detectable sequence specificity.

Gene ID:

ZNF346

Uniprot

Q9UL40

Synonyms:

Zinc finger protein 346 (Just another zinc finger protein), ZNF346, JAZ


Immunogen:

Recombinant Human Angiopoietin-related protein 4 protein (1-294AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: HeLa whole cell lysate, A549 whole cell lysate, MCF-7 whole cell lysate, 293T whole cell lysate, Jurkat whole cell lysate, PC-3 whole cell lysate. All lanes: ZNF346 antibody at 1:1000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 33, 36, 30 kDa. Observed band size: 33 kDa.