

ZNF346 Antibody, HRP conjugated



PACO64036

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds with low affinity to dsDNA and ssRNA, and with high affinity to dsRNA, with no detectable sequence specificity.

Gene ID:

ZNF346

Uniprot

Q9UL40

Synonyms:

Zinc finger protein 346 (Just another zinc finger protein), ZNF346, JAZ

Immunogen:

Recombinant Human Angiopoietin-related protein 4 protein (1-294AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A