

# ZNF346 Antibody, Biotin conjugated



PACO64038

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds with low affinity to dsDNA and ssRNA, and with high affinity to dsRNA, with no detectable sequence specificity.

**Gene ID:**

ZNF346

**Uniprot**

Q9UL40

**Synonyms:**

Zinc finger protein 346 (Just another zinc finger protein), ZNF346, JAZ

**Immunogen:**

Recombinant Human Angiopoietin-related protein 4 protein (1-294AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A