

IFRD1 Antibody, Biotin conjugated



PACO64050

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Could play a role in regulating gene activity in the proliferative and/or differentiative pathways induced by NGF. May be an autocrine factor that attenuates or amplifies the initial ligand-induced signal.

Gene ID:

IFRD1

Uniprot

O00458

Synonyms:

Interferon-related developmental regulator 1 (Nerve growth factor-inducible protein PC4), IFRD1

Immunogen:

Recombinant Human Interferon-related developmental regulator 1 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A