

RPIA Antibody, HRP conjugated



PACO64052

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, intracellular membrane-bounded organelle, identical protein binding, ribose-5-phosphate isomerase activity, D-ribose metabolic process, pentose-phosphate shunt, pentose-phosphate shunt, non-oxidative branch.

Gene ID:

RPIA

Uniprot

P49247

Synonyms:

Ribose-5-phosphate isomerase (EC 5.3.1.6) (Phosphoriboisomerase), RPIA, RPI

Immunogen:

Recombinant Human Ribose-5-phosphate isomerase protein (1-311AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A