

FABP4 Antibody, Biotin conjugated



PACO64074

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Lipid transport protein in adipocytes. Binds both long chain fatty acid, and retinoic acid, Delivers long-chain fatty acid, and retinoic acid, to their cognate receptors in the nucleus.

Gene ID:

FABP4

Uniprot

P15090

Synonyms:

Fatty acid, binding protein, adipocyte (Adipocyte lipid-binding protein) (ALBP) (Adipocyte-type fatty acid, binding protein) (A-FABP) (AFABP) (Fatty acid, binding protein 4), FABP4

Immunogen:

Recombinant Human Fatty acid, binding protein adipocyte protein (6-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A