

FGGY Antibody, Biotin conjugated



PACO64082

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

D-ribulokinase activity, neuron cellular homeostasis, pentose metabolic process.

Gene ID:

FGGY

Uniprot

Q96C11

Synonyms:

FGGY carbohydrate kinase domain-containing protein (EC 2.7.1), FGGY

Immunogen:

Recombinant Human FGGY carbohydrate kinase domain-containing protein (1-244AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A