

# HTR3E Antibody, Biotin conjugated



PACO64086

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. This receptor is a ligand-gated ion channel, which when activated causes fast, depolarizing responses. It is a cation-specific, but otherwise relatively nonselective, ion channel.

**Gene ID:**

HTR3E

**Uniprot**

A5X5Y0

**Synonyms:**

5-hydroxytryptamine receptor 3E (5-HT3-E) (5-HT3E) (Serotonin receptor 3E), HTR3E

**Immunogen:**

Recombinant Human 5-hydroxytryptamine receptor 3E protein (26-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A