

BHLHE23 Antibody, Biotin conjugated



PACO64166

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May function as transcriptional repressor. May modulate the expression of genes required for the differentiation and/or maintenance of pancreatic and neuronal cell types. May be important for rod bipolar cell maturation.

Gene ID:

BHLHE23

Uniprot

Q8NDY6

Synonyms:

Class E basic helix-loop-helix protein 23 (bHLHe23) (Class B basic helix-loop-helix protein 4) (bHLHb4), BHLHE23, BHLHB4

Immunogen:

Recombinant Human Class E basic helix-loop-helix protein 23 protein (1-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A