

# CALR Antibody, FITC conjugated



PACO64177

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Calcium-binding chaperone that promotes folding, oligomeric assembly and quality control in the endoplasmic reticulum (ER) via the calreticulin/calnexin cycle. This lectin interacts transiently with almost all of the monoglucosylated glycoproteins that are synthesized in the ER. Interacts with the DNA-binding domain of NR3C1 and mediates its nuclear export. Involved in maternal gene expression regulation. May participate in oocyte maturation via the regulation of calcium homeostasis.

**Gene ID:**

CALR

**Uniprot**

P27797

**Synonyms:**

Calreticulin (CRP55) (Calregulin) (Endoplasmic reticulum resident protein 60) (Erp60) (HACBP) (grp60), CALR, CRTC

**Immunogen:**

Peptide sequence from Human Calreticulin protein (185-196AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A