

CLDN11 Antibody, Biotin conjugated



PACO64206

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

Gene ID:

CLDN11

Uniprot

O75508

Synonyms:

Claudin-11 (Oligodendrocyte-specific protein), CLDN11, OSP OTM

Immunogen:

Recombinant Human Claudin-11 protein (23-82AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A