

Hmgb1 Antibody



PACO64271

Product Information

Size:

50ul

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:**Gene ID:**

Hmgb1

Uniprot

P63159

Synonyms:

High mobility group protein B1 (Amphoterin) (Heparin-binding protein p30) (High mobility group protein 1) (HMG-1), Hmgb1, Hmg-1 Hmg1


Immunogen:

Peptide sequence from Rat Hmgb1 protein (2-18AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: Rat Brain tissue. All lanes: Hmgb1 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 25 kDa. Observed band size: 30 kDa.