

Mmp2 Antibody



PACO64303

Product Information

Size:

50ul

Reactivity:

Mouse, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

Gene ID:

Mmp2

Uniprot

P33434

Synonyms:

72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) [Cleaved into: PEX], Mmp2


Immunogen:

Recombinant Mouse 72 kDa type IV collagenase protein (412-662AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: HeLa whole cell lysate. All lanes: Mmp2 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 75, 66 kDa. Observed band size: 75 kDa.