

Mmp2 Antibody, Biotin conjugated



PACO64306

Product Information

Size:

50ul

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Mmp2

Uniprot

P33434

Synonyms:

72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) [Cleaved into: PEX], Mmp2

Immunogen:

Recombinant Mouse 72 kDa type IV collagenase protein (412-662AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A