

# OR10A2 Antibody, Biotin conjugated



PACO64358

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Odorant receptor.

**Gene ID:**

OR10A2

**Uniprot**

Q9H208

**Synonyms:**

Olfactory receptor 10A2 (HP4) (Olfactory receptor OR11-86), OR10A2, OR10A2P

**Immunogen:**

Peptide sequence from Human Olfactory receptor 10A2 protein (247-262AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A