

OR5K2 Antibody



PACO64407

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

Odorant receptor.

Gene ID:

OR5K2

Uniprot

Q8NHB8

Synonyms:

Olfactory receptor 5K2 (Olfactory receptor OR3-9), OR5K2


Immunogen:

Peptide sequence from Human Olfactory receptor 5K2 protein (6-24AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: U251 whole cell lysate, U87 whole cell lysate. All lanes: OR5K2 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 36 kDa. Observed band size: 36 kDa.