

OR5V1 Antibody



PACO64419

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

Odorant receptor.

Gene ID:

OR5V1

Uniprot

Q9UGF6

Synonyms:

Olfactory receptor 5V1 (Hs6M1-21) (Olfactory receptor OR6-26), OR5V1

Immunogen:

Peptide sequence from Human Olfactory receptor 5V1 protein (302-318AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: A549 whole cell lysate. All lanes: OR5V1 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 37 kDa. Observed band size: 37 kDa.