

OR5V1 Antibody, Biotin conjugated



PACO64422

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Odorant receptor.

Gene ID:

OR5V1

Uniprot

Q9UGF6

Synonyms:

Olfactory receptor 5V1 (Hs6M1-21) (Olfactory receptor OR6-26), OR5V1

Immunogen:

Peptide sequence from Human Olfactory receptor 5V1 protein (302-318AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A