

PLPPR1 Antibody, Biotin conjugated



PACO64474

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Integral component of plasma membrane, nucleoplasm, lipid phosphatase activity, phosphatidate phosphatase activity, phospholipid dephosphorylation, phospholipid metabolic process, signal transduction.

Gene ID:

PLPPR1

Uniprot

Q8TBJ4

Synonyms:

Phospholipid phosphatase-related protein type 1 (Lipid phosphate phosphatase-related protein type 1) (Plasticity-related gene 3 protein) (PRG-3), PLPPR1, LPPR1 PRG3

Immunogen:

Peptide sequence from Human Phospholipid phosphatase-related protein type 1 protein (51-66AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A